



MERRY CHRISTMAS & HAPPY NEW YEAR

CLEARANCE SALE - GLOVES

Only while stocks last!

SafeClean Go Green Biodegradable Nitrile Gloves

A cleaner, greener and sustainable solution for a better environment.

Enviro-friendly | Biodegradable

The Go Green gloves are tested as per internationally recognised ASTM D5511 to determine their biodegradability in an anaerobic biodegradation environment.



\$92.⁹⁵

1 Carton (10xBox/100)

CODE	SIZE	CODE	SIZE
138100A	X-Small	138100C	Medium
138100B	Small	138100D	Large

Isolation Gown - SoftMed

- Level 1
- Soft Cuff
- Velcro neck, waist ties
- Non woven fabric
- Blue
- Made in Australia



Softmed

\$99.⁰⁰

Box/60



CODE	SIZE	CODE	SIZE
SM-SG301S	Small	SM-SG301XL	X-Large
SM-SG301M	Medium	SM-SG301XXL	XX-Large
SM-SG301L	Large		

Black Nitrile Ultra Soft Examination Gloves - Bastion

- Powder Free
- Nitrile, non latex
- Black colour
- Textured fingertips



\$114.²⁵

1 Ctn (10 x Box/100)

CODE	SIZE
BNG7853	Medium

Suretrax Process Indicator Labels - Meditrax



GETINGE
GETINGE GROUP

\$39.⁴⁰

Roll/700

CODE	COLOUR
PIPRP902A/BL	Blue

Latex Lightly Powdered Gloves - Bastion

- Beaded Cuff
- Excellent strength and tactility
- Cornstarch powder used is free of GE material
- Silicone free
- Micro Textured finish allows for outstanding wet & dry grip




\$125.²⁵

1 Ctn (10 x Box/100)

CODE	SIZE	CODE	SIZE
BNG2522	Small	BNG2524	Large
BNG2523	Medium		

Surgery Log Book - Meditrax



\$38.⁵⁰

Each

GETINGE
GETINGE GROUP

CODE	SIZE
SSMOO/6	N/A

Complete Peridontal / Examination Kit - MDDI

Each disposable sterile kit contains:

- 1 x Front Face Mirror
- 1 x Perio Probe/Sickle Double End
- 1 x Tweezers
- 2 x Face Masks with Ear Loop, Level 2
- 1 x Cotton Roll #2
- 1 x Tissue Wipe
- 1 x Tray



\$77.⁹⁰

Box/25 Kits



CODE	SIZE
SI550ST	N/A

Titems Chemical Indicators

Class 6

(Time, temperature, steam & pressure)

Changes colour from purple to green when all critical sterilisation parameters have been achieved



\$69.⁷⁵

Pkt/250

GETINGE
GETINGE GROUP

CODE	SIZE
AUS508	N/A

MASKS

Level 2

Masks - MGuard

- Earloop, Blue
- High filtration 3 layer mask.
- Conforms to AS4381:2015
- Deep pleats hold the mask off the face for improved breathing comfort.
- PFE 0.1 micron >98%
- BFE 3 micron >98%



\$7.⁹⁵
Box/50

CODE 70204

Level 3

Masks - Softskin Safe+ Medicom

- Blue, Tie-On
- Conforms to AS4381:2015
- Fluid Resistant 160mm Hg
- Softness and superior comfort
- Soft tissue inner layer to reduce skin irritation
- Anti-fog non-woven barrier on inner layer to prevent fogging



\$10.³⁵
Box/50



CODE SFS02050BF
STYLE Anti-Fog

PARTICULATE SURGICAL RESPIRATOR P2/N95

P2/N95 A-MED Particulate Surgical Respirator with headstraps - SoftMed

- P2, 95%
- Level 3, 160mmHg
- Fish style folded
- Elastic 2 x headloops
- Nose strip with foam
- Latex Free



\$49.⁹⁵
Box/20



CODE SM-RA201/20
STYLE Headloop

Particulate Respirator P2/N95 Flat Folded Earloop - SoftMed

- Pkt/10 Individually wrapped
- Latex Free
- Fluid Resistance 160mmHg



\$18.⁹⁵
Box/10



CODE SM-SR201/10
STYLE Earloop

S7XTRA HOSPITAL GRADE DISINFECTANT

ALCOHOL-FREE. KILLS COVID-19.
MANUFACTURED IN AUSTRALIA.

S-7XTRA: Multi-Purpose Biocidal Cleaning-Disinfectant Solution

- S-7XTRA is a multi-purpose biocidal cleaner that is also a hospital-grade 2-in-1 cleaner disinfectant.
- It cleans practically any surface or material, eradicating pathogens.
- With market-leading efficacy, S-7XTRA kills bacteria, viruses, fungi, spores, mycobacteria (such as TB) and moulds up to 99.9999% - in clean and dirty conditions.
- S-7XTRA has Reactive Barrier Technology (RBT) which means it remains active between cleaning cycles for ongoing protection.
- It has TGA-certified claims for a 24-hour residual activity on touched surfaces for SARS-CoV-2 (COVID-19), bacteria, and yeast. It also has certified claims for a 72-hour biostatic activity against SARS-CoV-2 and bacteria.
- S-7XTRA solutions come in a concentrate and hospital-grade Ready To Use (RTU) format giving both flexibility and certainty.
- Steri-7 also have a range of RTU impregnated wipes in either a pillow pack format (suited to hard-surfaces) or a tub format (suitable for equipment and where moisture needs to be minimised).
- S-7XTRA cleaner-disinfectants are listed with claims on the ARTG, including against SARS-Cov-2.



ARTG 232011 ARTG 232014

Wipes

Tub/200 Wipes
(wipes also available in pillow packs)

\$18.⁷⁵
Pkt/200

CODE S7TW-200
STYLE Wipes

Solution

S7XTRA Trigger Spray 750ml

\$14.⁹⁵
Each

CODE S7RTU-750ML
STYLE Spray

Solution

S7XTRA Ready To Use 5 Litre

\$39.⁸⁰
Each

CODE S7RTU-5LTR
STYLE RTU 5L

Solution

S7XTRA Concentrate 5 Litre
2 Pumps will dilute to make 500ml
3 Pumps will dilute to make 750ml

\$208.⁶⁰
Each

CODE S7CON-5LTR
STYLE Con 5L

ZHERMACK

Elite HD+ Impression Material

Each packet contains: 2 x 50ml cartridges + 6 mixing tips

\$63.³⁵

Per pkt

Tray Material / Heavy Body

CODE 90200



Regular Body

CODE 90205



Light Body Fast Set

CODE 90206



Light Body Normal Set

CODE 90207



Bite Registration

Occlufast + Colour - Zhermack

- Colour Change: Green to Yellow/Lime
- Detail Reproduction: 5 µm
- Working Time: 30 seconds
- Setting Time: 0:60 (min:s).
- High dimensional stability, up to 15 days
- Imperceptible consistency
- High final hardness (37±5 Shore D)
- Thixotropic

\$74.³⁵

2 x 50ml

Zhermack
BETTER INNOVATION

CODE 90282

Mouthguard / Denture Boxes - Large, Vented

**\$1.⁶⁵**Each
MIN ORDER 10

CODE	COLOUR	CODE	COLOUR
85041	Red	85045	Purple
85042	Green	85046	Black
85043	Hot Pink	85047	Yellow
85044	Orange	85049	Mixed

Alginates

Hydrogum 5 - Long Life Alginate

- Highly accurate 5 day dimensional stability alginate, easy to mix.
- Extra fast set -45 sec in mouth. Ideal for gaggers.
- Purple colour
- Berry flavoured, dust free

**\$14.⁰⁰**

1 x 453g

Zhermack
BETTER INNOVATION

CODE 90100

Orthoprint - Long Life Alginate

- Extra rapid set alginate, 45 second set in the mouth
- Ideal for orthodontics
- Yellow colour
- Pleasant vanilla flavor
- Easy to mix
- Highly elastic, dust free

**\$14.⁰⁰**

1 x 500g

Zhermack
BETTER INNOVATION

CODE 90116

Tropicalgin - Long Life Alginate

- Chromatic:
- Red during mixing
- Orange during the working phase
- Yellow during positioning in the patients mouth
- Mango Flavour
- Fast Set, 1 min set in mouth, dust free

**\$14.⁰⁰**

1 x 453g

Zhermack
BETTER INNOVATION

CODE 90105

Impression Trays - Disposable

Impression Trays

- Perforated
- Multi-Coloured Translucent with Handle
- Can be heat altered

**\$35.⁴⁵**

Per Pkt

Impression Trays
Perfo Orange SU

CODE	QTY
71201	48/Pkt

Impression Trays
Perfo Orange SL

CODE	QTY
71202	48/Pkt

Impression Trays
Perfo Blue MU

CODE	QTY
71203	48/Pkt

Impression Trays
Perfo Blue ML

CODE	QTY
71204	48/Pkt

Impression Trays
Perfo Yellow LU

CODE	QTY
71205	40/Pkt

Impression Trays
Perfo Yellow LL

CODE	QTY
71206	40/Pkt

Pink Upper Anterior

CODE	QTY
P1	48/Pkt

Pink Lower Anterior

CODE	QTY
P2	48/Pkt