

## SAFETY DATA SHEET

### Section 1: IDENTIFICATION

#### BACTEX CONCENTRATE

**Recommended Use:** Hospital Grade Disinfectant  
**Product Code:** 020206 (2x5L)

**Whiteley Corporation Pty. Ltd.**

(A.C.N. 000 906 678)

Postal Address: P. O. Box 1076 North Sydney NSW 2059

Telephone Number: (02) 9929 9155 Facsimile: (02) 9929 9077

Web: [www.whiteley.com.au](http://www.whiteley.com.au)

Emergency Telephone Number: Poisons Information Centre (National) 131126

Distributed in New Zealand by:

**Whiteley Corporation NZ Limited**

Address: PO Box 22-519 Otahuhu, Auckland, New Zealand 1640

Telephone Number: 092737313, 0800 257 352

Emergency Telephone Number: New Zealand Poisons Centre 0800 POISON or 0800 764 766

### Section 2: HAZARDS

#### GHS Classification

Skin Irritation (Category 2)  
 Eye Corrosion (Category 1)  
 Acute Toxicity - Oral (Category 4)

Signal Word        **DANGER**

#### Hazard Statements

H302	Harmful if swallowed
H315	Causes skin irritation
H318	Causes serious eye damage

#### Precautionary Statements

P260	Do not breathe dusts or mists
P264	Wash hands thoroughly after using
P270	Do not eat, drink or smoke when using this product.
P280	Wear protective gloves/protective clothing and eye protection

#### Response Statements

P301 + P312 + P330 + P331	IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. Call a POISON CENTRE or doctor/physician if you feel unwell.
P303 + P361 + P353	IF ON SKIN: Take off immediately all contaminated clothing. Rinse skin with water.
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310	Immediately call a POISON CENTRE or doctor/physician.
P363	Wash contaminated clothing before reuse

Disposal Statements

P501	Dispose of container as per local regulations
------	---

HSNO Classifications

6.1D	Substances that are acutely toxic – Harmful
6.3A	Substances that are irritating to the skin
8.3A	Substances that are corrosive to ocular tissue

**Section 3: COMPOSITION INFORMATION**

Ingredient	CAS No	Proportion
Benzalkonium Chloride	63449-41-2	10-30%
Ingredients deemed not to be hazardous	Not applicable	To 100%

**Section 4: FIRST AID**

<b>Eye (Contact)</b>	Hold eyelids apart and flush the eye continuously with running water. Immediately call a Poison Centre or doctor/physician.
<b>Skin (Contact)</b>	Remove contaminated clothing and flush skin and hair with running water. If skin irritation occurs seek medical advice. Wash contaminated clothing before reuse.
<b>Inhalation (Breathing)</b>	Remove to fresh air.
<b>Ingestion (Swallowing)</b>	DO NOT induce vomiting. Give water to drink if conscious. Contact a Poisons Information Centre (Phone 131126) or a doctor/physician.
<b>Advice to Doctor</b>	Treat symptomatically.
<b>First Aid Facilities</b>	Ensure an eye wash is available and ready for use.
<b>Additional Information</b>	No aggravated medical conditions are known to be caused by exposure to this product.

**Section 5: FIREFIGHTING MEASURE**

<b>Suitable Extinguishing Media</b>	Solution does not burn. Use extinguishing media suited to the materials that are burning. E.g. Dry chemical, CO <sub>2</sub> or water spray.
<b>Hazards from Combustion Products</b>	Carbon dioxide, carbon monoxide and other toxic gases may be produced in the case of fire or during thermal decomposition.
<b>Precautions for Fire Fighters and Special Protective Equipment</b>	Fire-fighters must wear full protective clothing including self-contained breathing apparatus and chemical splash suit. Ensure that no spillage enters drains or water courses. Remove from the vicinity containers not involved in the fire.
<b>Additional Information</b>	<b>Hazchem Code</b> –N/A

**Section 6: ACCIDENTAL RELEASE MEASURES**

<b>Emergency Procedure</b>	SAA/SNZ HB76: Dangerous Goods – Initial Emergency Response Guide (Guide 37) – Not Applicable.
----------------------------	---

**Spills / Clean up**

Clean up personnel should wear full protective clothing. Restrict access until completion of clean up. Then ensure adequate ventilation. Stop leak if safe to do so. Contain spill with absorbent material, such as towelling, sand, vermiculite or other inert material. Prevent spill entering stormwater drains or waterways. Collect and dispose of clean up material according to local regulations. Wash away remnants with copious amounts of cold water to sewer. Clean area by working from the periphery to the centre of spill or from the edge of the room to the centre.

<b>Section 7: HANDLING AND STORAGE</b>
--

**Precautions for Safe Handling**

Contact Whiteley Corporation sales representative for advice when using this product for any application other than that outlined on the label or technical bulletin.

Any non-intended or non-authorized use of this product may result in personal injury or damage to equipment.

Store product in original container.

Wash thoroughly after handling product.

**Conditions for Safe Storage**

Store in a cool, dry, well ventilated area. Keep container tightly sealed. Store below 30°C.

<b>Section 8: EXPOSURE CONTROL/PERSONAL PROTECTION</b>
--

**National Exposure Standards** – Source: National Exposure Standards for Atmospheric Contaminants in the Occupational Environment [NOHSC:1003] & New Zealand Workplace Exposure Standards 2002.

Ingredient

None known

CAS NoES-TWAES-STEL**Biological Limit Values**

No data available.

**Engineering Controls**

Use only in a well ventilated area.

**Personal Protective Equipment**

Eye/face protection – Safety glasses / face shield / chemical resistant goggles should be worn to prevent eye contact.

Skin protection – Use Nitrile gloves or similar to prevent skin contact.

Respiratory protection – Respirator is not usually necessary but if required use a half face filter respirator suitable for organic vapours.

Section 9: PHYSICAL AND CHEMICAL PROPERTIES			
---	--	--	--

<b>Appearance</b>	Slightly thick, clear blue liquid	<b>Odour</b>	Characteristic
<b>Odour Threshold</b>	No data available	<b>pH</b>	7.50 – 8.50
<b>Freezing Point</b>	Approx. 0°C	<b>Initial Boiling Point</b>	Approx. 100°C
<b>Boiling Range</b>	No data available	<b>Flash point</b>	No data available
<b>Evaporation Rate</b>	No data available	<b>Flammability</b>	No data available
<b>Upper flammability limit</b>	No data available	<b>Lower flammability limit</b>	No data available
<b>Vapour Pressure</b>	No data available	<b>Vapour Density</b>	No data available
<b>Relative Density</b>	1.03	<b>Solubility</b>	Completely miscible with water
<b>Partition Coefficient: n-octanol/water</b>	No data available	<b>Autoignition temperature</b>	No data available
<b>Decomposition temperature</b>	No data available	<b>Viscosity</b>	Approx. 100 cPs

Section 10: STABILITY AND REACTIVITY	
--------------------------------------	--

<b>Chemical Stability</b>	Product is stable and will not undergo any hazardous reactions under normal conditions of use and storage.
<b>Conditions to avoid</b>	None known.
<b>Incompatible materials</b>	None known
<b>Hazardous decomposition products</b>	None known.
<b>Hazardous reactions</b>	None known.

Section 11: TOXICOLOGICAL INFORMATION	
---------------------------------------	--

<b>Acute toxicity</b>	Benzalkonium Chloride: LD <sub>50</sub> 240mg/kg (Oral, rat)
<b>Skin corrosion/irritation</b>	Product causes irritation, pain and reddening on skin contact. Serious burns may result if the affected area if product is not removed immediately by thorough washing with water.
<b>Serious eye damage/irritation</b>	Product causes irritation, pain and reddening on eye contact. Serious, permanent eye damage may result if not treated immediately.
<b>Respiratory or skin sensitisation</b>	No data available.

<b>Germ Cell Mutagenicity</b>	No data available.
<b>Carcinogenicity</b>	No data available.
<b>Reproductive toxicity</b>	No data available.
<b>Specific Target Organ Toxicity – Single Exposure</b>	No data available.
<b>Specific Target Organ Toxicity – Repeated Exposure</b>	No data available.
<b>Aspiration Hazard</b>	No data available.

<b>Section 12: ECOLOGICAL INFORMATION</b>
---

<b>Ecotoxicity</b>	No data available.
<b>Persistence and degradability</b>	No data available.
<b>Bioaccumulative potential</b>	No data available.
<b>Mobility in soil</b>	No data available.
<b>Other adverse effects</b>	No data available.

<b>Section 13: DISPOSAL CONSIDERATIONS</b>
--

<b>Disposal method</b>	Disposal to sewer is normally recommended with copious amounts of water. Refer to State/Territory Land Waste Management Authorities if applicable. Containers are recyclable and can be disposed of by a licensed waste contractor. Containers can be disposed of to general waste or rinsed thoroughly and recycled.
<b>Special precautions</b>	Suitable for incineration by approved agent.

<b>Section 14: TRANSPORT INFORMATION</b>
--

Not classified as Dangerous Goods by the criteria of the Australian Dangerous Goods Code (ADG Code).	
<b>UN Number</b>	Not applicable.
<b>UN Proper Shipping Name</b>	Not applicable.
<b>Class and subsidiary risk</b>	Not applicable.
<b>Packing Group</b>	Not applicable.
<b>Special precautions for Hazchem Code</b>	Not applicable.

<b>Section 15: REGULATORY INFORMATION</b>
---

Poisons Schedule (SUSMP): schedule 6 – POISON.

All ingredients are listed in the Australia Inventory of Chemical Substances (AICS).

HSNO Approval Code: HSR002526

This document has been produced in accordance with the requirements of the Globally Harmonised System of Classification and Labelling.

<b>Section 16: OTHER INFORMATION</b>
--------------------------------------

**Date of preparation:** 22 April 2021

© This document remains property of Whiteley Corporation Pty Ltd.  
Alterations are not permitted without prior written authorisation from Whiteley Corporation Pty Ltd.