

OPAL[®] Instrument Grade - High Level Disinfectant



ARTG No: 206575

OPAL[®] is an Ortho-Phthalaldehyde (OPA) based high level disinfectant that features a reuse period of up to twice the length of competitors.

OPAL[®] is suitable for manual disinfection of compatible semi critical medical and dental devices as well as automated reprocessing of endoscopes and accessories.

OPAL[®] has the same contact time and temperature requirements for Manual and Automated Reprocessing. Because OPAL[®] is a room temperature high level disinfectant at 20 °C, there is no heating of the product

necessary for any reprocessing of devices. The 6-minute contact time with OPAL[®] for both Manual and Automated Reprocessing means faster turnaround times, saving time and money.

Packsizes available:

Size	Code
6 x 1L	150106
2 x 5L	150103

For more information on OPAL high level disinfectant or to register for a product trial and inservicing in your facility contact marketing@whiteley.com.au or call our Product Support Hotline on **1800 833 566**.