

# V-Wipes™

INSTRUMENT GRADE DISINFECTANT (LOW LEVEL)

ARTG: 354053

**V-Wipes™ are an Instrument Grade Disinfectant (low level) wipe, which means they are suitable for the disinfection of environmental and non-critical medical devices.**

The V-Wipes™ unique formula with its high quality apertured fabric is proven to kill a variety of bacteria and viruses including SARS-CoV-2 [COVID-19], *Enterococcus faecalis* [VRE]; *Klebsiella pneumoniae* [CPE/CRE], Measles virus and *Staphylococcus aureus* [MRSA].

V-Wipes™ unique pH neutral formula allows for fast, effective cleaning and disinfection of all surfaces within the healthcare environment. V-Wipes™ are intended for use on environmental surfaces and for non-critical medical devices such as hospital beds and theatre trolleys.

**V-Wipes™ Benefits:**

- ✓ Class IIb Instrument Grade Disinfectant - Low Level Wipes.
- ✓ pH neutral – safe for use on most surfaces & equipment.
- ✓ Highly absorbent apertured fabric.
- ✓ Leaves surfaces streak free.
- ✓ Extremely durable.
- ✓ Non hazardous.
- ✓ Does not contain phenolics or chlorine.
- ✓ Non flammable.
- ✓ Friendly to use on metals [non corrosive].
- ✓ PHMB free



**Product codes:**

Flatpack – 12x80 Wipes (230x330mm): 210634  
 Canister – 4x100 Wipes (180x250mm): 210561



**Learn more or find a distributor**

<b>Australia</b>	ABN No 24 000 906 678	<b>New Zealand</b>	NZBN 9429041224243
 <b>1800 833 566</b>		 <b>0800 257 352</b>	
 <b>www.whiteley.com.au</b>		 <b>www.whiteley.co.nz</b>	
PO Box 1076, North Sydney NSW 2059		PO Box 22-519, Otahuhu, Auckland 1640	