



# NEW PRODUCT

## FINGERTIP PULSE OXIMETER

A small device intended for measuring the pulse oxygen saturation, pulse rate and perfusion index through a users finger.

Suitable for monitoring of adult and paediatric patients in medical clinics, hospitals and at home.

### Features

- OLED display of SpO2, PR, PI, pulse bar and waveform
- Suitable for adult and paediatric patients
- Special splash proof and drop resistant design
- Spot check and continuous measurement mode
- Pulse rate analysis for spot check measurement
- Up to 12 groups SpO2 data storage
- Audible & visual alarm with low battery indication
- Automatic power on/off
- Four direction display
- Artefact removal and anti-motion
- Accurate and sensitive



### Specifications:

<b>SpO2</b>	
Transducer	Dual-wavelength LED sensor
Measuring range	35-100%
Measuring accuracy	≤ 2% range from 70-100%
<b>Pulse Rate</b>	
Measuring range	30-240bpm
Measuring accuracy	±2bpm or ±2% (whichever is the greater)
<b>Perfusion Index Display Range</b>	
0-20%	
<b>Default Alarm Limit</b>	
SpO2 low limit	90%
Pulse rate high limit	120bpm
Pulse rate low limit	50bpm
<b>Measuring Modes</b>	
30 seconds spot check mode Continuous mode	
<b>Size</b>	56 x 34 x 30mm (L x W x H)
<b>Weight</b>	52g
<b>Power</b>	2 x AAA alkaline battery

Order Code: PMFPO

### Pacific Medical

Phone: 1300 076 332  
Email: [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)  
[www.pacificmedical.com.au](http://www.pacificmedical.com.au)

